



# LikoStretch Mod 600 IC



## For Extra-Demanding Situations

This stretcher is specifically developed for horizontal lifting and transferring of patients with very high demands on careful handling. This can be, for instance, severely burnt patients or patients with serious fractures in the neck or back who consequently cannot be turned to the side.

The weight-bearing lift straps are applied gently using a spatula, without having to lift or roll the patient. LikoStretch 600 IC can also be used together with a full LiftSheet for simpler transfer situations.



## ■ ACCESSORIES



### Stretch Leveller

Recommended accessory enabling easy adjustment of the centre of gravity when the patient is lifted in LikoStretch 600 IC.

Fits: 3156065 and 3156065B

### Mast Attachment Prod. No. 2001600 (not illustrated)

Mast Attachment Golvo for LikoStretch 600 IC.

4 straps incl.

Recommended when LikoStretch 600 IC is used in combination with mobile lift Golvo in situations requiring increased stability during the lifting motion.

## ■ PRODUCT OVERVIEW

3156065	LikoStretch Mod 600 IC. (250 kg/550 lbs) Standard 70 cm (28 inch) excl. accessories
3156065B	LikoStretch Mod 600 IC. (250 kg/550 lbs) Width 80 cm (32 inch) excl. accessories
3156066	LikoStretch Mod 600 IC. (200 kg/440 lbs) Adjustable, width 60-90 cm (24-36 inch) excl. accessories
3156200	Stretch Leveller
3166506	SideRails Mod 600 IC (pair)
3684001	LiftStrap, 8 pcs (set)
3684002	Spatula
3682504	LiftSheet 800, standard: 80 cm (32 inch), reinforced fabric, excl SideRails
3682508	LiftSheet 1800, width: 90 cm (36 inch), reinforced fabric, excl SideRails
3682509	LiftSheet 1800, standard: 80 cm (32 inch), reinforced fabric, excl SideRails
3166508	SideRails 800 (pair)
3166518	SideRails 1800 (pair)



LikoStretch 600 IC complies with the requirements for Medical Device Directive Class I products.

Maximum load: 200-250 kg (440-550 lbs)

Design and Quality by Liko in Sweden.

Liko is quality certified in accordance with ISO 9001 and its equivalent for the medical device industry, ISO 13485.

Liko is also certified in accordance with environmental standard ISO 14001.



[www.liko.com](http://www.liko.com)

#### Manufactured by:

Liko AB  
Nedre vägen 100  
975 92 Luleå  
Sweden  
[info@liko.se](mailto:info@liko.se)

#### Local representative in UK

Hill-Rom Limited  
Clinitron House, Ashby Park  
Ashby de la Zouch,  
Leicester LE65 1JG  
Tel: 01530 411000  
Fax: 01530 411555  
[info@liko.co.uk](mailto:info@liko.co.uk)

#### To find your local Liko representative,

please visit [www.liko.com](http://www.liko.com) or call  
+46 (0)920 47 47 00



A Hill-Rom Company