



# FlexoStretch



## Full flexibility, also for heavy patients

The FlexoStretch is intended for lifting and transfer of patients in a fully supine position. The body support part consists of either a lift sheet or five separate lift straps. The stretcher width can be adjusted to between 75 and 90 cm (29.5 and 35.5 inches) by moving the lift arms in or out.

FlexoStretch is equipped with a Stretch leveller that enables adjustment of the centre of gravity during lifting. Horizontal stretchers can be used together with a variety of Liko lifts.

Read more on our website: [www.liko.com](http://www.liko.com).



FlexoStretch is a medical device class 1 product.  
Max. load: 300 kg, 660 lbs.

Design and Quality by Liko in Sweden.

Liko is quality certified in accordance with ISO 9001 and its equivalent for medical devices, ISO 13485.  
Liko is certified in accordance with the Environmental Management System ISO 14001.



## PRODUCT DATA

3156057	FlexoStretch
3684100	LiftStrap XL (set of 5 straps; 2 pcs XL80 + 3 pcs XL100)
3684101	LiftStrap XL80
3684102	LiftStrap XL100
3684105	LiftSheet XL regular
3684106	LiftSheet XL tall

[www.liko.com](http://www.liko.com)

### Manufactured by:

Liko AB  
Nedre vägen 100  
975 92 Luleå  
Sweden  
[info@liko.se](mailto:info@liko.se)

### Local representative in UK

Hill-Rom Limited  
Clinitron House, Ashby Park  
Ashby de la Zouch,  
Leicester LE65 1JG  
Tel: 01530 411000  
Fax: 01530 411555  
[info@liko.co.uk](mailto:info@liko.co.uk)

To find your local Liko representative,  
please visit [www.liko.com](http://www.liko.com) or call  
+46 (0)920 47 47 00



A Hill-Rom Company