



UltraStretch



For heavy supine lifting

UltraStretch is intended for lifting and transfer of heavy patients in a fully supine position. UltraStretch is used in combination with Liko's lift system FreeSpan UltraTwin, with rails 80 cm (31.5 inches) apart, or with UltraTwin in a ceiling mounted installation.

UltraStretch is easily mounted to the lift units via its Quick-release Hooks. No tools are required. The body support part consists of either a lift sheet or five separate lift straps.

■ PROD. NO.

3156058	UltraStretch
3684100	LiftStrap XL (set of 5 straps; 2 pcs XL80 + 3 pcs XL100)
3684105	LiftSheet XL regular
3684106	LiftSheet XL tall



Maximum load: 400 kg, 880 lbs.

Design and Quality by Liko in Sweden.

Liko is quality certified according to ISO 9001 and its equivalence for the medical device industry, ISO 13485. Liko is also certified according to environmental standard ISO 14001.



www.liko.com

Manufactured by:

Liko AB
Nedre vägen 100
975 92 Luleå
Sweden
info@liko.se

Local representative in UK

Hill-Rom Limited
Clinitron House, Ashby Park
Ashby de la Zouch,
Leicester LE65 1JG
Tel: 01530 411000
Fax: 01530 411555
info@liko.co.uk

To find your local Liko representative,
please visit www.liko.com or call
+46 (0)920 47 47 00



A Hill-Rom Company