


Slings for all needs

BASIC MODELS FOR SAFE, EASY LIFTING




UNIVERSALSLING, MOD. 00/02

- Good basic model, ideal for most patients
- Provides support for the entire back
- The patient holds his/her arms inside the sling
- Enables an upright sitting posture (00)
- Mod. 02 allows extra lifting height

SIZE: S, MS, M, L, XL

FABRIC: Polyester, Net polyester, plastic coated net,




SOFT ORIGINAL HB. SLING, MOD. 26

- Can be used to lift most patients
- Best choice for lifting from a horizontal position
- Optimized fit and comfort
- Provides added comfortable support for the head
- Enables a semi-reclined sitting posture

SIZE: XS, S, MS, M, L, XL, XXL

FABRIC: Polyester, teddy bear patterned polyester jersey




ORIGINAL HIGHBACK SLING, MOD. 20/21/25

- Can be used to lift most patients
- Good choice for lifting from a horizontal position
- Provides support for the entire back and head
- Enables a semi-reclined sitting posture
- Mod. 25, Solo HighBack, is a disposable sling intended for short term use (do not wash!)

SIZE: XS, S, MS, M, L, XL, XXL

FABRIC: Net polyester, plastic coated net, polyester, teddy bear patterned jersey, nonwoven material




ORIGINALSLING, MOD. 10/11

- Allows for active patient participation
- Good for lifting from seated/semi-seated position
- Supports the back up to the shoulder blades
- The patient holds his/her arms inside or outside the sling
- Easy to apply

SIZE: XS, S, MS, M, L, XL

FABRIC: Polyester, net polyester, plastic coated net, teddy bear patterned jersey

MODELS FOR SPECIAL NEEDS




COMFORTSLING PLUS, MOD 300/310/350/360

- Good for leg amputees and patients who suffer from painful joints
- The sling is left under the patient
- Enables an open angle at the hips, in combination with 4-point slingbar
- Provides added comfortable support for the head (350, 360)

SIZE: M, L, XL

FABRIC: Net polyester, plastic coated net




AMPUTEESLING, MOD. 70/75

- For patients with single or double femoral (thigh) amputation
- Ideal for patients who tend to slide out of slings
- Provides support for the entire back and head (75)
- Semi-reclined sitting posture

SIZE: S, M, L, XL

FABRIC: Polyester, Net polyester




SILHOUETTESLING MOD. 22

- A pliable sling in net fabric for lifting to shells and moulded seats. Good for patients with disabilities
- Can be left in the chair after lifting
- Provides support for the entire back and head
- The patient holds his/her arms inside or outside the sling.

SIZE: XS, S, M, L

FABRIC: Net polyester




ULTRASLING, MOD. 06

- Designed for bariatric patients
- Allows plenty of extra space for the pannus
- Can be used for toileting
- Max. load 500 kg/1,100 lbs

SIZE: Short, Regular, Long

FABRIC: Polyester, (Net polyester)

STAND/AMBULATE/MOBILITY




MASTERVERST, MOD. 60

- For patients who can stand on their own
- Safe and secure for raising, standing, walking and balance training
- Allows natural movement during raising
- Good for toileting, lower body free accessible
- Can be used to advantage with Leg Harness

SIZE: XXS, XS, S, M, L, XL, XXL

FABRIC: Polyester, plastic coated net, teddy bear patterned polyester jersey




MASTERVERST, MOD. 64

- Same application area as Mod. 60
- Features front fastening
- Front-adjustable for bust and stomach
- For independent patients
- Leg Harness (accessory) may be suitable for some patients

SIZE: XS, S, M, L, XL

FABRIC: Polyester, plastic coated net, teddy bear patterned polyester jersey




LIFTPANTS, MOD. 92/920*

- For patients who can stand on their own
- Safe and secure for raising, standing, ambulate and balance training
- Allows natural movement during raising
- Ultra LiftPants, Mod. 920 developed for ambulate training of bariatric patients

SIZE: XS, S, M, L, XL, *XXL, *XXXL

FABRIC: Net polyester, teddy bear patterned polyester jersey




VEST FOR STANDINGSHELL, MOD. 67

- Lifts the patient in a standing shell from a horizontal to a standing position
- Valuable aid to habilitation of children and youths
- Leg Harness included
- Does not strain the standing shell

SIZE: XS, S, M, L

FABRIC: Polyester, teddy bear patterned polyester jersey

TOILET/HYGIENE




HYGIENESLING, MOD. 40/41

- For patients with good stability
- Easy to apply
- Enables dressing/undressing during lifting to/from the toilet
- The patient holds his/her arms outside the sling
- Mod. 41 is equipped with a belt for added safety and security

SIZE: XS, S, M, L, XL

FABRIC: Polyester, teddy bear patterned polyester jersey




HYGIENESLING, MOD. 45/46

- Equipped with Safety Belt, ideal for patients with poor stability
- Enables dressing/undressing during lifting to/from the toilet
- Can be used with one arm held inside the sling
- Mod. 46 features a back support for extra support, safety and comfort

SIZE: XS, S, M, L, XL

FABRIC: Polyester, teddy bear patterned polyester jersey




HYGIENEVEST, MOD. 50

- Appropriate for patients with high/low tonus and spasticity
- Enables dressing/undressing during lifting to/from the toilet
- Supports and lifts around the upper body
- High degree of safety when lifting to/from the toilet
- Upright sitting posture

SIZE: XXS, XS, S, M, L, XL, XXL, XXXL

FABRIC: Polyester, teddy bear patterned polyester jersey



HYGIENEVEST, MOD. 55

- Ideal for patients with spasticity
- Enables dressing/undressing during lifting to/from the toilet
- Provides support for the entire back and head
- High degree of safety when lifting to/from the toilet
- Supports and lifts around the upper body

SIZE: XXS, XS, S, M, L

FABRIC: Polyester, teddy bear patterned polyester jersey