



Solo Disposable Products



Products for short-term use

The spread of multi-resistant bacteria (MRB) within healthcare, along with the increased demand for products for short-term use, has led to the development of Solo – a range of disposable lifting accessories.

What is unique about most of the Solo products, compared to our ordinary lifting accessories, is that they are made from a very strong non-woven material. Solo RepoSheet is the exception, as this product is made in a cotton/polyester fabric. Solo products are not intended to be washed but are to be disposed of once they have become soiled or when the patient no longer needs them. In other words, they are designed for short-term use by one particular patient.

Apart from being an effective solution for preventing the spread of infection, the lifting accessories also have economical and practical advantages: laundering is eliminated and they are always immediately available on the ward.



1 Solo HighBack Sling™

Solo HighBack Sling is a sling that suits most patients. Solo HighBack Sling can be used in all lifting situations, for instance, between bed and wheelchair and to the toilet.

2 Solo RepoSheet™

Solo RepoSheet is used to move up patients who have slid down in the bed or to turn bed-ridden patients, for instance, in order to wash them or to avoid bedsores.

3 Solo Octo LiftSheet™

Solo Octo LiftSheet is intended for horizontal lifting with Liko's OctoStretch. The OctoStretch is easy to use and appreciated by surgical and intensive care departments.

4 Solo LiftSheet 1800™

Solo LiftSheet 1800 is intended for horizontal lifting with LikoStretch Mod 600 IC. LikoStretch Mod 600 IC can also be used together with washable lift straps when extra high demands are put on the lifting situation. The lift straps are placed under the patient without the need for turning the patient to the side.

5 Solo MultiStrap™

Solo MultiStraps are used to lift and hold parts of the patient's body, in order to facilitate examination, care and treatment. Solo MultiStraps can also be used to turn a patient in bed and to support the patient when turned to the side.

6 Solo SupportVest™

Solo SupportVest is used with the Sabina Sit-to-Stand Lift for raising a patient from a sitting to a standing position. Solo SupportVest lifts behind the patient's back and under the arms. It suits patients who need a lot of help to raise themselves but who can then stand firmly on their own legs.

7 Solo Cover for Ultra LiftPants™

Solo Cover is a disposable accessory for Ultra LiftPants. The cover is slipped over the LiftPants to protect them from becoming soiled.

Contact your Liko representative or visit www.liko.com for more information on our Solo products.



Design and Quality by Liko in Sweden

Solo disposable products are class 1 medical devices.

Liko is quality certified in accordance with ISO 9001 and its equivalent for the medical device industry, ISO 13485.

Liko is certified in accordance with environmental standard ISO 14001.



www.liko.com

Manufactured by:

Liko AB
Nedre vägen 100
975 92 Luleå
Sweden
info@liko.se

Local representative in UK

Hill-Rom Limited
Clinitron House, Ashby Park
Ashby de la Zouch,
Leicester LE65 1JG
Tel: 01530 411000
Fax: 01530 411555
info@liko.co.uk

To find your local Liko representative,
please visit www.liko.com or call
+46 (0)920 47 47 00



A Hill-Rom Company